

ALDEFREE

Instant Surface Disinfectant

Manufactured in GMP unit

Made in India

Product of
Indigeneous
Research
& Inputs

Other Range

Infectrol® GF

Ultrabasil® Plus

Enzycide®

Disdril

Remorust

Microbat®



UPS Hygienes Pvt. Ltd.

(AN ISO 9001-2015 Certified Company)
Customer Care No.: 022-65214567 / 9821338361
E-mail: info@upshygiene.com

Your Partner In Infection Control

Aldehyde Free, Alcoholic Disinfectant Pre Injection Spray

- OT** - On Operation table, Lights, Monitors, For emergency disinfection in between surgeries.
- ICU** - ECG machines, Monitors, Trolleys, Ventilator.
- AKD** - Respirators, Monitors, Haemodialysis units
- NICU** - Incubators, Ventilators, Monitors, Pumps, Cables, Dental Chair, CT Scanners, Pharmaceutical sterile sections.

For preparing dressing, skin disinfection, before and after injection for infants & diabetics, before dressing minor injuries & sutures

Composition

Benzalkonium Chloride 0.3 w/v
2-propanol 70% v/v

Packing

100 ml	60 bott./carton
250ml + Spray	30 bott./carton
5 ltrs.	4 cans/carton

BKC goes into deeper layer giving long term microbiocidal effect.
Alcohol degreases the skin.

- Alcohol 70% has excellent germicidal activity against gram+ve & gram-ve bacteria including MRSA, VRE, M.tuberculosis & fungi^{1,2}.
- Alcohol 70% is effective against ROTAVIRUS³.
- Benzalkonium Chloride 0.1-0.3% is effective against most of the bacteria^{4,5}. It is antifungal.

EFFECTIVE AGAINST BACTERIA, FUNGI, VIRUSES LIKE HIV & ROTA.

Dries fast within 10-25 seconds.

Do not switch on or off the electronic devices while spraying

Ref.:

1. Journal of Infectious diseases 1940, 67, 275-81
2. Pharmaceutical Journal 1936, 137, 79-81.
3. Applied Environmental microbiology, 1989, 55, 3113-8
4. Principle & Practises of Disinfection by Russel Hugo.
5. Preservation & Sterilisation by Sykes G. 3rd edition

Website : www.upshygiene.co.in

Helpline : 022-65214567 / 9821338361

