



UPS Hygienics Pvt. Ltd.

(AN ISO 9001-2015 Certified Company)

Customer Care No.: 022-65214567

E-mail: info@upshygiene.com

ULTRABASIL[®] 5Q

Multipurpose Disinfectant Concentrate

Composition :-

Alkyl Dimethyl benzyl ammonium chloride	8.19%
Didecyl dimethyl ammonium chloride	8.70%

**5th Generation
Quarternary
Ammonium
Compound.**

Bactericidal, Fungicidal, Virucidal, Algacidal.

- Gram + ve organisms- S. aureus, E. hirae, L.monocytogenes.
- Gram -ve organisms-P.aeuroginosa, E. coli, S. typhimurium, C. jejuni.
- Fungus- C. albicans, A. niger.

※ One-Step Disinfectant Cleaner

Dilution:

- For critical areas: 4-8 ml in 1 litre of water (critical area)
- High risk areas like OT, ICU, ICCU, NICU, and AKD.
- For non critical areas: 2 ml in 1 litre of water
- Areas like reception, wards, path labs etc.

ENVIRONMENT / SURFACES / FLOORING / INSTRUMENT